Massage

JOHNSON®

I

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Thank you for choosing a **JOHNSON** product. **Massager** is designed to provide a deep muscle massage, promote relaxation, and improve muscle recovery.

Please read this manual carefully before use and keep it for future reference.

GB

GENERAL WARNINGS

WARNING: Consult your doctor before use if: you have an open wound, are taking blood thinners or have other serious vascular problems, have heart problems or use a pacemaker or defibrillator, are pregnant or have a serious medical condition, or want to massage areas of the body with minimal muscle coverage.

Do not use this device to massage sensitive areas such as the head, face, cervical region, or vertebral/spinal cord.

WARNING: This appliance can be used by children aged 8 years and above and persons with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

Cleaning and user maintenance shall not be made by children unless they are older than 8 years and supervised. Keep the appliance and its cord out of reach of children under 8 years of age.

Children shall not play with the appliance.

Failure to follow these instructions or improper use of the product may cause damage or injury.

This appliance, suitable for non-professional use only, must be used only for the purpose for which it was specifically designed, as explained in this manual. Any other use is considered improper and therefore dangerous, and will void the warranty.

The manufacturer cannot be held responsible for any damage resulting from improper, incorrect, or unreasonable use. Installation must be carried out according to the manufacturer's instructions. Incorrect installation may cause damage to people, animals, or property; the manufacturer cannot be held responsible.

Do not leave the appliance exposed to atmospheric agents (frost, rain, sun, etc.). Do not keep the appliance near heat sources (e.g., radiators).

Packaging materials (plastic bags, nails, ex-

panded polystyrene, etc.) must be kept out of the reach of children or infirm persons, as they represent potential sources of danger. We also remind you that these must be disposed of separately.

Before using the device for the first time, make sure you remove all protective labels or sheets except those containing technical data or product information.

If the device has air vents, make sure they are never blocked, even partially, and that no objects such as needles, pins, or similar are inserted, and that water does not enter.

The device should not be used after a fall or if it is damaged.

Use the device only with the accessories provided. To avoid compromising safety, never use a damaged accessory; have it checked, repaired, or replaced with an approved one.

Never immerse the device or its cord in water or any other liquid. Do not rinse under the tap. Use a damp cloth for cleaning.

Do not insert foreign objects into the device. Keep cables away from traffic areas or where there is a risk of tripping or falling.

Do not place objects on top of the device.

Always use the original USB cable provided in

the box.

Do not touch the cable or the device connected to the USB port with wet or damp hands.

Always disconnect the device from the USB port before cleaning, maintenance, or moving the device.

The device must be powered only by a safety extra-low voltage corresponding to the marking on the product.

Before use, check the USB cable for visible damage. If the cable is damaged, replace it with an approved one. Turn off the device before inserting the USB cable to charge.

Be sure to unplug the USB cable when the battery is fully charged.

We recommend using it in dry environments. If using it in a humid environment, make sure that water never comes into contact with the device.

Do not place the device near open flames or heat sources, and do not use it in dusty environments or in the presence of flammable vapors. The device contains a rechargeable battery: once charged, the device can be used without a power connection.

Do not place the USB cable under rugs, cloths, or similar objects.

It is strictly forbidden to use the device while it is charging. Recharge the battery using the provided USB cable. Be sure to unplug it when charging is complete.

Never touch moving parts.

Do not remove components, screws, etc., and do not attempt to disassemble or repair the device. The device automatically turns off after a 15 minute operating cycle, and it is not recommended to use the massager on the same muscle area again.

This product is not a medical device and is not intended for the diagnosis, treatment, or prevention of diseases and/or pathologies. During use, keep your beard, hair, etc., away from the moving parts of the device to avoid injury.

HOW TO DISPOSE OF THE DEVICE

When disposing of this device, the battery must be removed and it is also recommended that it be rendered inoperable.

WARNING: Batteries must be disposed of safely; they must be placed in appropriate containers or handed over to the environmental recycling department for disposal. It is also recommended that any potentially hazardous

parts of the device be rendered harmless.

WARNINGS ON USING RECHARGEABLE BATTERIES

If battery fluid comes into contact with skin or eyes, rinse the affected area with water and consult a doctor.

Risk of ingestion: Children can swallow rechargeable batteries and suffocate. Therefore, keep rechargeable batteries out of reach of children.

Before inserting batteries, observe the positive (+) and negative (-) polarities.

If battery fluid leaks, wear protective gloves and clean the battery compartment with a dry cloth. Protect rechargeable batteries from excessive heat and cold.

Risk of explosion: Do not throw rechargeable batteries into a fire. Do not disassemble, open, or crush rechargeable batteries.

Fully charge rechargeable batteries before using them for the first time.

To achieve the longest possible battery life, fully charge rechargeable batteries at least twice a year.

WARNING:

When the battery reaches the end of its life (it can no longer be recharged), the device must be disposed of as waste.

WARNING

If the product contains glass parts, they are not covered under the warranty. Plastic parts are not covered under the warranty.

WARNING

Damage to the power cord (if present) resulting from wear and tear is not covered under the warranty; any repairs will be the responsibility of the owner.

WARNING

If it becomes necessary to take or ship the appliance to an authorized service center, we recommend thoroughly cleaning all parts. If the appliance is even slightly dirty or encrusted, or has limescale, dust, or other deposits, externally or internally, for health and hygiene reasons, the service center will reject the appliance without examining it.

The symbol  on the product or packaging indicates that the product must not be treated as household waste, but must be taken to the appropriate collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise arise from inappropriate waste handling of this product. For more detailed information about recycling this product, please contact your local city office, your household waste disposal service, or the shop where you purchased the product.

The  symbol on the product or on its packaging indicates that this battery may not be treated as household waste. Instead, it should be taken to the applicable collection point for the recycling of used batteries.

For more detailed information about recycling this battery, please contact your local city office, your household waste disposal service, or the shop where you purchased the product.

COMPONENTS

- A. On/off button and speed adjustment
- B. USB-C connector for charging cable
- C. LED charging and operation lights
- D. Massage head attachment
- E. Spherical head
- F. U-shaped head
- G. Flat head
- H. Cylindrical head
- L. USB-C charging cable
- M. Head housing

Remove the appliance parts from the packaging and dispose of the packaging according to legal requirements. Check that the product is intact and that all parts are present. If there are obvious anomalies or damage during transport, do not use the appliance; contact your retailer or an authorized service center.

The percussion massager is used to relieve muscle tension and temporarily increase circulation, reducing the risk of injury. Percussion pressure also helps eliminate lactic acid buildup from muscles and reduce pain.

USE

Before first use, fully charge the device, wash and thoroughly dry all parts that will come into contact with the body.

Charging

Connect the USB-C cable (L) to the USB-C charging connector (B) and then to a USB charging source.

During charging, the LEDs will flash blue. When charging is complete, the LEDs will stop flashing and remain solid blue. Once charging is complete, disconnect the USB-C cable (L) from the source.

A full charge takes approximately 4 hours, depending on the power of the USB charging source.

A full charge allows up to approximately 20 use cycles of 15 minutes each (battery life varies depending on the intensity selected and the pressure applied during use).

WARNING: It is recommended to keep the device fully charged to always have maxi-

mum power available.

Note: If the battery charge drops below 5% during use, the LEDs will begin flashing, depending on the selected speed. At this point, the device needs to be recharged.

To check the battery charge level before using the percussion massager, quickly press the on/off and speed control button (A). The LEDs will light up according to the detected battery charge level (1 = almost empty, 5 = fully charged).

WARNING: Battery life and performance can be affected by various factors, including usage time, environmental conditions, and charging and discharging cycles. Over time and with prolonged use, battery performance may naturally decrease.

Massage Head Assembly

Select the massage head based on the area to be treated and your personal needs. The device comes with four different interchangeable heads:

Spherical head

for treating large muscle groups, such as the buttocks, legs, and arms, and larger areas of the body, such as the pectorals and lower back.

U-shaped head

for areas near the spine and the shoulder blades. Suitable for relaxing and shaping the entire body.

Cylindrical head

ideal for specific areas, such as the palms of the hands and feet, and for targeting the deep nerve structures in the joints.

Flat head

for a broad and extensive massage, suitable for all areas of the body.

Once you have selected the accessory you want to use, insert the head coupling (D) into the housing (M) on the device, applying light pressure. The head coupling (D) must be fully inserted into the housing (M).

To remove the head, pull it outward until it comes completely loose.

Tip:

When removing the accessory, it is best to pull it and turn it counterclockwise, as if it were a threaded connection, to make it easier to detach.

When attaching the accessory, it is best to push it and turn it clockwise, as if it were a threaded connection, to make it easier to attach.

Caution: Make sure the device is turned off before inserting or removing a massage head.

Use

Always turn the device on before placing it on the area to be treated.

By holding down the on/off and speed control button (A) for a couple of seconds, The device turns on and operates at the lowest massage intensity (1).

By repeatedly pressing the on/off and speed control button (A), you can increase the massage intensity from level 1 (lowest) to level 5 (highest). The blue LED indicators (C) above the on/off and speed control button (A) light up to indicate the selected intensity (e.g., massage intensity 3 LEDs lit in blue: 3).

To turn off the device, move it away from the treated area and press and hold the on/off and speed control button (A) again for a couple of seconds.

The device turns off automatically after 15 minutes of continuous use to avoid excessive muscle stimulation.

To resume treatment, repeat the start-up sequence. It is advisable to start treatment with the lowest massage intensity and then gradually increase it if necessary to avoid possible injury or trauma. Move the device gently, applying light pressure to the area to be treated, taking care not to linger on the same spot for more than 60 seconds. Use the device on clean, dry skin, free of clothing and other objects that could hinder the movement of the massage head.

N.B.: The device cannot be turned on while charging.

HOW TO USE A PERCUSSION MASSAGER

Always start with the largest muscle groups, such as the glutes, quadriceps, calves, and back, until you have mastered the technique correctly and safely.

Use your other hand as a guide to identify the bone sections to avoid.

Be careful not to press too hard or stay in one spot for too long. If used incorrectly, the percussion massager can easily cause more harm than good. Remember to always keep the percussion massager moving for best results.

Use the percussion massager for short periods, focusing on each part of the body for a few minutes, to get the maximum benefit from the treatment without causing trauma or discomfort to the area.

Be sure to avoid using the percussion massager on arteries, nerves, and internal organs;

the area above the neck, armpits, and other areas with less muscle coverage, as well as areas near organs and the aorta.

On bones and joints, and on any bruised or injured body part.

The percussion massager should only be used on soft muscle tissue without causing pain or discomfort.

Pay attention to your body's response, and if you feel pain or discomfort, stop treatment immediately.

CLEANING AND MAINTENANCE

Make sure the device is turned off and the USB-C cable is disconnected from the charging source before proceeding with cleaning and maintenance.

Use a slightly damp cloth to clean the device and dry with a soft, dry cloth.

Always avoid allowing moisture to enter the device's case.

When storing or moving, keep it in the provided case.

Note: If the percussion massager is not used for an extended period of time, we recommend charging it at least once a month. This will keep the battery in optimal condition.

A discharged battery will reduce its lifespan and compromise performance.

CAUTION: Do not immerse the device in water. Keep away from liquids and/or heat sources.

WARNING:

When the battery reaches the end of its life (it can no longer be recharged), the device must be disposed of as waste.

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GUARANTEE

The validity of the guarantee starts from the purchasing date certified by the seller's stamp and expires at the end of the indicated period, even when the device has not been used. The guarantee covers all substitutions or repairs due to material or manufacturing defects. The guarantee has no validity in case of device tampering or when the defect is due to incorrect or improper use of the same. Substitutions and repair covered by the guarantee are made free of any charge for goods delivered free our service company representative. Only the freight will be at buyer's expense. All replaced parts or devices will become our property.

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MOD. Massage

DATA ACQUISTO
PURCHASE DATE

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TIMBRO CENTRO ASSISTENZA
AFTER SALES SERVICE CENTER STAMP

MASS-0426-I

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JOHNSON® S.r.l